



Thermometry
1-800-962-9888



COVIDIEN

positive results for life™

Table of Contents

Thermometry

TYMPANIC THERMOMETERS

Genius™ 2 Infrared Tympanic Electronic Thermometer	4
--	---

ELECTRONIC THERMOMETERS

Filac™ 3000 AD Electronic Thermometer	5
Filac™ 3000 EZ Electronic Thermometer.....	6
Filac™ FasTemp™ Electronic Thermometer	7

TYMPANIC THERMOMETER PROBE COVERS

Genius™ 2 Tympanic Electronic Thermometer Probe Covers....	8
Genius™ Model 3000A and FirstTemp™ 3000 Tympanic Thermometer Probe Covers.....	8

ELECTRONIC THERMOMETER PROBE COVERS

Filac™ 3000 Electronic Thermometer Probe Covers	8
Filac™ FasTemp™ Electronic Thermometer Probe Covers.....	9

ACCESSORIES

Mobile Thermometer Stands	9
Calibration Units	9
Locking Bases.....	10
Side and Top Utility Baskets.....	11
Calibration Plugs	11
Isolation Chambers	11

Tympanic Thermometers



Genius™ 2 Infrared Tympanic Thermometer

- 1-2 second response
- Tapered probe design
- Peak Select System analyses up to 100 separate measurements
- No touch disposable probe cover application and ejection
- Push button alternates Fahrenheit/Celsius
- 60 second pulse timer
- Takes 15,000 temperature reading using 3 AAA batteries
- On base probe cover storage
- Last temperature recall
- ASTM & CEN compliant
- Utilizes item 303030 disposable probe cover

CODE	DESCRIPTION	SHIP CASE
303000	Genius 2 Thermometer with Base	1
303099	Genius 2 Calibration Checker	1

Genius 2 Thermometer hardware is warranted for one year from shipment for manufacturer's defect only.

Genius 2 Calibration Checker hardware is warranted for one year from shipment for manufacturer's defect only.

Electronic Thermometers



Filac™ 3000 AD Electronic Thermometer

- Quick Mode: Accuracy: +/- 0.3°C, Speed: 3-4 seconds
- Standard Mode: Accuracy: +/- 0.1°C, Speed: 6-10 seconds
- Direct Mode: Accuracy: +/- 0.1°C, Speed: <60 seconds
- Cold Mode for hypothermic patients
- Icon based user interface on backlit LCD display
- Electronic anti theft system
- Push button alternates Fahrenheit/Celsius
- Mountable base with key lock
- Interchangeable, color coded isolation chambers with matching probes
- Measurement range 30°C to 43°C (86°F to 109°F)
- Takes 6000 temperature reading using 4 AA batteries
- "Auto-On" "Auto-off" feature
- Last temperature recall
- 15, 30, 45, and 60 second pulse timer
- Utilizes item 202020, 502000, or 500500 disposable probe cover

CODE	DESCRIPTION	SHIP CASE
505000	Filac 3000 AD Oral/Axillary Complete	1
505005	Filac 3000 AD Rectal Complete	1
505027	Filac 3000 AD Oral/Axillary Complete - 9ft	1
505037	Filac 3000 AD Rectal Complete - 9ft	1
500026	Filac 3000 Oral Probe - 4ft	1
500027	Filac 3000 Oral Probe - 9ft	1
500028	Filac 3000 Oral Isolation Chamber	1
500036	Filac 3000 Rectal Probe - 4ft	1
500037	Filac 3000 Rectal Probe - 9ft	1
500038	Filac 3000 Rectal Isolation Chamber	1
500099	Filac 3000 Calibration Plug	1

Filac 3000 AD hardware is warranted for one year from shipment for manufacturer's defect only.



Filac™ 3000 EZ Electronic Thermometer

- Quick Mode: Accuracy: +/- 0.3°C, Speed: 3-4 seconds
- Standard Mode: Accuracy: +/- 0.1°C, Speed: 6-10 seconds
- Direct Mode: Accuracy: +/- 0.1°C, Speed: <60 seconds
- Icon based user interface on backlit LCD display
- Push button alternates Fahrenheit/Celsius
- Interchangeable, color coded isolation chambers with matching probes
- Measurement range 30°C to 43°C (86°F to 109°F)
- Takes 6000 temperature reading using 4 AA batteries
- "Auto-On" "Auto-off" feature
- Last Temperature recall
- Utilizes item 202020, 502000, or 500500 disposable probe cover

CODE	DESCRIPTION	SHIP CASE
504000	Filac 3000 EZ Oral/Axillary Complete	1
504005	Filac 3000 EZ Rectal Complete	1
504027	Filac 3000 EZ Oral/Axillary Complete - 9ft	1
504037	Filac 3000 Rectal Complete - 9ft	1
500026	Filac 3000 Oral Probe - 4ft	1
500027	Filac 3000 Oral Probe - 9ft	1
500028	Filac 3000 Oral Isolation Chamber	1
500036	Filac 3000 Rectal Probe - 4ft	1
500037	Filac 3000 Rectal Probe - 9ft	1
500038	Filac 3000 Rectal Isolation Chamber	1
500099	Filac 3000 Calibration Plug	1

Filac 3000 EZ hardware is warranted for one year from shipment for manufacturer's defect only.

Electronic Thermometers



Filac™ FasTemp™ Electronic Thermometer

- 4 - 10 second reading (oral, axillary, rectal)
- Color coded isolation chambers with matching probes
- Global iconography
- Oversized backlit LCD display
- Push button alternates Fahrenheit/Celsius
- Last temperature recall
- 30 second pulse timer
- Temperature range (85°F - 110°F/30°C - 43°C)
- Electronic anti-theft system
- One year warranty
- ASTM and CEN compliance
- 9 volt alkaline battery included
- Utilizes item 202020, 502000, or 500500 disposable probe cover

CODE	DESCRIPTION	SHIP CASE
202000	Complete Oral/Axillary System, including Electronic Thermometer, 4 Language Operating Manual Blue Oral/Axillary Isolation Chamber, Blue Probe with 4ft Cord, Wall Mountable Base, Mounting Hardware	1
202018	FasTemp Monitor Only	1
202026	Blue Oral/Axillary Probe with 4ft Cord	1
202028	Blue Oral/Axillary Isolation Chamber	1
202034	Red Rectal Probe with 4ft Cord	1
202036	Red Rectal Isolation Chamber	1
202038	Red Rectal Isolation Chamber, Red Rectal Probe with 4ft Cord, Wall Mountable Base, Mounting Hardware	1
202027	Blue Oral/Axillary Probe with 9ft Cord	1
202037	Red Rectal Probe with 9ft Cord	1
202099	FasTemp Calibration Plug	1

Filac FasTemp hardware is warranted for one year from shipment for manufacturer's defect only.

Tympanic Thermometer Probe Covers



Genius™ 2 Tympanic Thermometer Probe Covers

- Disposable, single-use
- No touch probe cover application and ejection
- Easy load probe cover
- Smooth, rigid plastic with clear membrane
- Two probe cover cassettes fit into base of thermometers

CODE	DESCRIPTION	SHIP CASE
303030	Genius 2 Probe Covers Cassette of 16 Probe Covers 6 cassettes/package; 11 packages/box; 2 boxes per case	2,112 each/case



Genius™ Model 3000A Tympanic Thermometer Probe Covers

- Disposable, single-use
- Smooth, rigid plastic with clear membrane
- No touch probe cover application and ejection
- Probe cover cassette fits into base of thermometers
- May also be used by FirstTemp™ 2000A tympanic thermometer

CODE	DESCRIPTION	SHIP CASE
8884810055	Cassette of 21 Probe Covers; 5 cassettes/package; 10 packages/box; 2 boxes per case	2,100 each/case

Electronic Thermometer Probe Covers



Filac™ 3000 AD and EZ Electronic Thermometer Probe Covers

- Disposable, single-use
- Smooth, rigid plastic
- Probe cover package fits into isolation chamber
- Ejected into waste receptacle by depressing release button
- May also be used by Filac FasTemp™ electronic thermometer

CODE	DESCRIPTION	SHIP CASE
500500	Filac 3000 FastTemp Probe Covers 20 probe covers; 25 packages/box; 1 box per case	500 each/case
502000	Filac 3000 FastTemp Probe Covers 20 probe covers; 25 packages/box; 1 box per case	2,000 each/case

Electronic Thermometer Probe Covers



Filac™ FasTemp™ Electronic Thermometer Probe Covers

- Disposable, single-use
- Smooth, rigid plastic
- Probe cover package fits into isolation chamber
- Ejected into waste receptacle by depressing release button
- May also be used by Filac™ 3000 AD and EZ electronic thermometers

CODE	DESCRIPTION	SHIP CASE
202020	Filac FasTemp Probe Covers 20 probe covers; 25 packages/box; 10 boxes per case	5,000 each/case

Tympanic Thermometer Featured Accessories



Mobile Stand

- Portable thermometer stand with probe cover holder
- Includes Genius™ 2 tympanic thermometer locking mount
- Compatible with the Filac™ 3000 mounting base for Mobile Stand, item 500047

CODE	DESCRIPTION	SHIP CASE
303059	Genius 2 Mobile Stand with locking mount	1



Genius™ 2 Cart Conversion Kit

- Bracket to attach Genius 2 thermometer to rolling stand models 8884830301 and 202057

CODE	DESCRIPTION	SHIP CASE
303060	Genius 2 Cart Conversion Kit	1



Genius™ 2 Calibration Checker

- Checks calibration status of Genius tympanic thermometer

CODE	DESCRIPTION	SHIP CASE
303099	Genius 2 Calibration Checker	1

Infrared Tympanic Thermometer Featured Accessories



Genius™ 2 Locking Wall Mount

- Mobile thermometer stand with side utility basket
- Helps reduce risk of theft

CODE	DESCRIPTION	SHIP CASE
303058	Genius Locking Wall Mount	1

Filac™ 3000 Electronic Thermometer Featured Accessories



Filac™ 3000 Mounting Base with Pole Clamp

- Locking base with pole clamp for attachment to mobile stands

CODE	DESCRIPTION	SHIP CASE
500048	Filac 3000 Mounting Base with Pole Clamp	1



Filac™ 3000 Mounting Base with Lock

- Locking base to secure Filac 3000 thermometer to the wall

CODE	DESCRIPTION	SHIP CASE
500046	Filac 3000 Wall Mounting Base with Lock	1



Filac™ 3000 Mounting Base for Rolling Stand

- Mounting base to secure Filac 3000 thermometer to mobile stand
- Filac 3000 mounting base with lock not included

CODE	DESCRIPTION	SHIP CASE
500047	Filac 3000 Mounting Base for Mobile Stand	1

Filac™ 3000 Electronic Thermometer Featured Accessories



Filac™ 3000 Calibration Plug

- Calibration plug to check calibration of the Filac 3000 thermometer


CODE	DESCRIPTION	SHIP CASE
500099	Filac 3000 Calibration Plug	1



Utilities Basket for Mobile Stand

- Basket with pole clamp for mobile stand

CODE	DESCRIPTION	SHIP CASE
8884813711	Thermometer Mobile Stand Utility Basket	1



Filac™ 3000 Thermometer Probes

- Replacement probes for Filac 3000 thermometer
- Available in Oral/Axillary or Rectal versions
- Available in 4ft or 9ft cord lengths

CODE	DESCRIPTION	SHIP CASE
500026	Filac 3000 Oral Probe 4ft	1
500027	Filac 3000 Oral Probe 9ft	1
500036	Filac 3000 Rectal Probe 4ft	1
500037	Filac 3000 Rectal Probe 9ft	1



Filac™ 3000 Isolation Chamber

- Available in Oral/Axillary or Rectal versions
- Isolation chamber in house thermometer probes
- Does not include Filac 3000 thermometer probe

CODE	DESCRIPTION	SHIP CASE
500028	Filac 3000 Oral Isolation Chamber	1
500038	Filac 3000 Rectal Isolation Chamber	1

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
202000	Complete Oral/Axillary System,	7	500028	Filac 3000 Oral Isolation Chamber	5, 6, 11
202018	FasTemp™ Monitor Only	7	500036	Filac 3000 Rectal Probe - 4ft	5, 6, 11
202020	Filac™ FasTemp™ Probe Covers	9	500037	Filac 3000 Rectal Probe - 9ft	5, 6, 11
202026	Blue Oral/Axillary Probe with 4ft Cord	7	500038	Filac 3000 Rectal Isolation Chamber	5, 6, 11
202027	Blue Oral/Axillary Probe with 9ft Cord	7	500046	Filac 3000 Wall Mounting Base with Lock	10
202028	Blue Oral/Axillary Isolation Chamber	7	500047	Filac 3000 Mounting Base for Mobile Stand	10
202034	Red Rectal Probe with 4ft Cord	7	500048	Filac 3000 Mounting Base with Pole Clamp	10
202036	Red Rectal Isolation Chamber	7	500099	Filac 3000 Calibration Plug	5, 6, 11
202037	Red Rectal Probe with 9ft Cord	7	500500	Filac 3000 FasTemp Probe Covers	8
202038	Red Rectal Isolation Chamber, Red Rectal Probe, 4ft Cord, Wall Mountable Base, Mounting Hardware	7	502000	Filac 3000 FasTemp Probe Covers	8
202099	FasTemp Calibration Plug	7	504000	Filac 3000 EZ Oral/Axillary Complete	6
303000	Genius 2™ Thermometer with Base	4	504005	Filac 3000 EZ Rectal Complete	6
303030	Genius 2 Probe Covers	8	504027	Filac 3000 EZ Oral/Axillary Complete - 9ft	6
303058	Genius Locking Wall Mount	10	504037	Filac 3000 Rectal Complete - 9ft	6
303059	Genius 2 Mobile Stand with the locking mount	9	505000	Filac 3000 AD Oral/Axillary Complete	5
303060	Genius 2 Cart Conversion Kit	9	505005	Filac 3000 AD Rectal Complete	5
303099	Genius 2 Calibration Checker	4, 9	505027	Filac 3000 AD Oral/Axillary Complete - 9ft	5
500026	Filac 3000 Oral Probe - 4ft	5, 6, 11	505037	Filac 3000 AD Rectal Complete - 9ft	5
500027	Filac 3000 Oral Probe - 9ft	5, 6, 11	8884810055	Genius Model 3000A Cassette of 21 Probe Covers	8
			8884813711	Thermometer Mobile Stand Utility Basket	11



COVIDIEN, COVIDIEN with logo and™ marked brands are trademarks of Covidien AG or an affiliate.
© 2008 Covidien AG or an affiliates All rights reserved.

H6200TH 10M 7/08

15 HAMPSHIRE STREET
MANSFIELD, MA
02048

800-962-9888
508-261-8000

WWW.COVIDIEN.COM